

MPC Series Couplings

For medical device, biotech and pharmaceutical applications

Economical media connections

Polycarbonat and polysulfon MPC couplings replace expensive glass and stainless steel fittings. The bio-compatible materials meet USP Class requirements and pass all physicochemical, cytotoxicity, leachables and extracebles testing. Using disposable MPC series quick disconnect couplings eliminates recertification and resterilisation. Laboratory procedures are simplified with fewer steps. GMP/GLP and FDA compliance is easier to maintain. It all adds up to cost savings.




STEDIM's Flexboy® media bay
pigtail connections use MPC series
couplings

Sterilize MPC couplings your way

Polysulfon is autoclavable to 25 cycles. Or use EtO, gamma radiation or E-beam methods. Easy to use push-button latch reduces handling of tubing and fittings that may contaminate the fluid path. The silicone O-ring seal always makes a reliable, leak free connection.

Bioprocessing applications

- ◆ Cell culture and reagent packaging
 - ◆ Sterile fluid
 - ◆ Bio-reactor lines
 - ◆ Transfer lines
 - ◆ Series connections
 - ◆ Delivery and recovery lines
 - ◆ Laboratory tubing connections
- ⇒ **New! Now single packed, sterile!**



	Sterilisation methods*				
	Autoclave Parts uncoupled	EtO	Gamma Radiation	E-Beam	
Polycarbonat (PC)	At 250 °F for 30 min.	-	-	-	* Colder Products Company did not conduct any sterilisation test above 5.0 megarads. Polycarbonate and Polysulfone manufacturer's material specifications may differ from actual MPC test results. EtO cycle parameters available upon request. Working pressures and other technical information have been prepared from actual test results and other data considered to be reliable.
	Yes, up to 10 cycles	Yes, up to 5 cycles	Yes*, up to 5.0 megarads	Yes*, up to 5 megarads	
	Polycarbonate is rated for 250 °F	-	Coupling color shift to clear at 2.5 megarads	Coupling color shift to faded purple tint	
Polysulfon (PS)	at 270 °F for 30 Min.	-	-	-	
	Yes, up to 25 cycles	Yes, up to 5 cycles	Yes*, up to 5.0 megarads	Yes*, up to 5.0 megarads	