



F\_01\_hp\_engl\_Biopure 1\_R01\_090409.doc




Warnings: Over 6meg rads some weakening and also discolouration may occur. Do not over torque – maximum torque is 2Nm. After autoclaving, Clamps may need to be re-torqued. Not to be used with solid PTFE gaskets / seals.




**Information**: These "Bio Barb Sanitary"-products are manufactured in an ISO Class 8 Clean Room. These products can be sterilised by gamma irradiation at a dose of 25 kGy or autoclave sterilised at temperature up to 134°C.



F\_01\_hp\_engl\_Biopure 2\_R01\_090409.doc




**Information**: These "Bio Barb Mini"-products are manufactured in an ISO Class 8 Clean Room. These products can be sterilised by gamma irradiation at a dose of 25 kGy or autoclave sterilised at temperature up to 134°C.




**Information**: These "Bio End Caps"-products are manufactured in an ISO Class 8 Clean Room. These products can be sterilised bay gamma irradiation at a dose of 25 kGy or autoclave serilised at temperature up to 134°C.



F 01 hp engl Biopure 3 R01 090409.doc



Warnings: Over 6meg rads some weakening and also discolouration may occur.



**Information**: These "Bio Y's"-products are manufactured in an ISO Class 8 Clean Room. These products can be sterilised by gamma irradiation at a dose of 25 kGy or autoclave sterilised at temperature up to 134°C.