We have a disposable clamp to be utilized on large bag tubing. The component is meant to seal off the fluid path on tubing with a simple unlocking feature to open the fluid path.

Typical applications that this large clamp might be used for consist of but are not limited to: biopharmaceutical manufacturing and pharmaceutical processes; disposable filtration systems; mixing containers and vessels; assemblies and tubing sets.

## Features/Benefits

Controls fluid flow

Variety of tubing sizes - from 3/4" to 1 1/2" (external tube dimensions) 1/8" wall


Secure lock

Reusable durability, disposable economy

Autoclave and irradiation stable\* Manufactured utilizing Makrolon 2458 WH Polycarbonate Resin

Ability to add a tamper evident tie

\*Recommendation – do not clamp tubing during irradiation or autoclave cycle, distortion may occur.



Preliminary Data Sheet Vers. 222060

**Applied Critical Fluids GmbH** Industriestr. 53, D-68169 Mannheim, Germany phone: 0049 (0)621/3 09 80 93-30 fax: 0049 (0)621/3 09 80 93-99 mobile: 0049/171/3572946 email: wz@acftec.de internet: www.acftec.de