

Applied Critical Fluids GmbH

Industriestrasse 53


68169 Mannheim

Tel.: +49 (0) 621 309 809 30

Fax.: +49 (0) 621 309 809 399

email: info@acftec.de

0



Features	Benefits
Ergonomic thumb latch	Easy to operate – even with gloved hands
USP Class VI materials	Meets biocompatibility requirements
Sterilizable by autoclave, Et0, e-beam, or gamma	Reusable, yet economical enough to allow disposability
Parting line-free hose barb	Eliminates potential leak path

disconnects. In addition, coupling halves can be rotated when

connected to reduce tube kinks.

Specifications

Pressure:

Vacuum to 60 psi, 4.14 bar

Temperature:

Polysulfone: -40° F to 300° F (-40° C to 149° C) Polycarbonate: -40° F to 250° F (-40° C to 121° C)

ABS: -40° F to 160° F (-40° C to 71° C)

Materials:

Main components:

Polycarbonate (purple tint), USP Class VI Polysulfone (amber tint), USP Class VI ABS (white), USP Class VI

Locking sleeves:

Polysulfone (white), (not applicable for ABS)

O-rings: Silicone (clear), platinum-cured, USP Class VI

and Buna-N (black), USP Class V

Sterilization:

Gamma: Up to 50 kGy irradiation


Polycarbonate: 250°F (121°C), 30 minutes, up to

10 repetitions. Sterilize uncoupled only.

Polysulfone: Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.

Tubing Sizes: 1/4" to 3/8" ID, 6.4mm to 9.5mm ID

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of Colder's products in their own application conditions.



This graph is intended to give you a general idea of the performance capabilities of each product line. The shaded area of the graph represents the operating range of the product family, i.e. upper and lower values are shown. Therefore, depending on the exact coupling configurations selected, you can reasonably expect values to fall within the shaded area

mates with SaniQuik Sanitary Series (See pages 8-9)

Liquid Flow Rates

Liquid Flow Rate Information for Couplings

The chart below shows the flow rate for Colder couplings. Each coupling was tested with water at 70° F (21° C). To determine flow rates for specific coupling configurations use the formula below.

$$Q = C_V \sqrt{\frac{\Delta P}{S}}$$


Q=Flow rate in gallons per minute C_v =Average constant of various rates (see chart) Δ P=Pressure drop across coupling (psi)

S=Specific gravity of liquid

C_v Values for MPC Couplings

BODIES	MPC 22004T	MPC 22006T
MPC170004T	2.8	2.8
MPC170006T	2.8	5.5

Product Dimensions



Coupling Bodies



ABS

TERMINATION	TUBING	METRIC EQ.	FLOW	STRAIGHT THRU	Α	В
IN-LINE	1/4" ID	6.4mm ID	.210"	MPC17004T	.93 (23.6)	1.30 (33.0)
HOSE BARB	3/8" ID	9.5mm ID	.290"	MPC17006T	.93 (23.6)	1.30 (33.0)



POLYCARBONATE

TERMINATION IN-LINE HOSE BARB	TUBING 1/4" ID 3/8" ID	METRIC EQ. 6.4mm ID 9.5mm ID	FLOW .210" .290"	STRAIGHT THRU MPC17004T03 MPC17006T03	A .93 (23.6) .93 (23.6)	B 1.30 (33.0) 1.30 (33.0)	
TERMINATION IN-LINE HOSE BARB WITH LOCK	TUBING 1/4" ID 3/8" ID	METRIC EQ. 6.4mm ID 9.5mm ID	FLOW .210" .290"	STRAIGHT THRU MPCK17004T03 MPCK17006T03	A .99 (25.2) .99 (25.2)	B 1.30 (33.0) 1.30 (33.0)	



POLYSULFONE

TERMINATION IN-LINE HOSE BARB	TUBING 1/4" ID 3/8" ID	METRIC EQ. 6.4mm ID 9.5mm ID	FLOW .210" .290"	STRAIGHT THRU MPC17004T39 MPC17006T39	A .93 (23.6) .93 (23.6)	B 1.30 (33.0) 1.30 (33.0)	
TERMINATION IN-LINE HOSE BARB WITH LOCK	TUBING 1/4" ID 3/8" ID	METRIC EQ. 6.4mm ID 9.5mm ID	FLOW .210" .290"	STRAIGHT THRU MPCK17004T39 MPCK17006T39	A .99 (25.2) .99 (25.2)	B 1.30 (33.0) 1.30 (33.0)	

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

Accessories



DESCRIPTION

Leash plug for MPC body

Leash cap for MPC insert

MATERIALS

Soft, flexible, medical-grade PVC PART NO.

MPC30L

MPC32L

Soft, flexible, medical-grade PVC Note: For validation quantities of MPC and MPX, contact Colder for 25 piece bag quantities

DID YOU KNOW ...

Colder's products for Life Sciences applications are manufactured in our ISO Class 7 certified cleanroom. The SMC, MPC, MPX, MPU, Sanitary, HFC39 and Steam-Thru® product lines are all molded from medical-grade materials and are packaged in double bags with material certifications.

DID YOU KNOW ...

Many of Colder's connectors are made from Animal-Free materials thereby reducing the amount of BSE-related documentation. Contact Customer Service at 1-800-444-2474 or 651-645-0091 for further information about Colder's Animal-Free material offering.

Coupling Inserts ABS

TERMINATION IN-LINE HOSE BARB	TUBING 1/4" ID 3/8" ID	6.4mm ID 9.5mm ID	.210" .290"	STRAIGHT THRU MPC22004TM MPC22006TM	O-RING Silicone Seal USP Class VI Silicone Seal USP Class VI	A .60 (15.2) .60 (15.2)	B 1.30 (33.0) 1.30 (33.0)
TERMINATION	TUBING	METRIC EQ.	FLOW	STRAIGHT THRU	O-RING	` ,	B
IN-LINE	1/4" ID	6.4mm ID	.210"	MPC22004T	Buna-N Seal USP Class V		1.30 (33.0)
HOSE BARB	3/8" ID	9.5mm ID	.290"	MPC22006T	Buna-N Seal USP Class V		1.30 (33.0)



POLYCARBONATE

TERMINATION	TUBING	METRIC EQ.	FLOW	STRAIGHT THRU	O-RING	` ,	B
IN-LINE	1/4" ID	6.4mm ID	.210"	MPC22004T03M	Silicone Seal USP Class VI		1.30 (33.0)
HOSE BARB	3/8" ID	9.5mm ID	.290"	MPC22006T03M	Silicone Seal USP Class VI		1.30 (33.0)
TERMINATION	TUBING	METRIC EQ.	FLOW	STRAIGHT THRU	O-RING	` ,	B
IN-LINE	1/4" ID	6.4mm ID	.210"	MPC22004T03	Buna-N Seal USP Class V		1.30 (33.0)
HOSE BARB	3/8" ID	9.5mm ID	.290"	MPC22006T03	Buna-N Seal USP Class V		1.30 (33.0)





POLYSULFONE

TERMINATION IN-LINE		METRIC EQ. 6.4mm ID	FLOW .210"	STRAIGHT THRU	O-RING Silicone Seal USP Class VI	A 60 (15.2)	B 1 30 (33 0)
HOSE BARB	-,	9.5mm ID	.290"		Silicone Seal USP Class VI	` '	` '



SEALING CAP	SEALING CAP W/LOCK	MATERIAL	A	B
MPC32003	MPCK32003	Polycarbonate	.93 (23.6)	1.30 (33.0)
MPC32039	MPCK32039	Polysulfone	.99 (25.2)	





SEALING PLUG	O-RING	MATERIAL	A	B
MPC30003M	Silicone Seal USP Class VI	Polycarbonate	.75 (19.1)	1.24 (31.5)
SEALING PLUG	O-RING	MATERIAL	A	B
MPC30039M	Silicone Seal USP Class VI	Polysulfone	.75 (19.1)	1.24 (31.5)

Accessories



All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

